



HQExoTM **Exosome-serum**

Catalog: Exo-BF04

PRODUCT INFORMATION

Name	HQExo™ Exosome-serum
Cat No.	Exo-BF04
Source	Exosome derived from human serum (small EVs) (healthy donors)
Product Overview	Exosomes are small membrane vesicles with a size of 30-160 nm that are released by different cell types. They can be found in various body fluids including plasma, malignant ascites, urine, amniotic fluid and saliva, which contains a tissue-specific signature wherein a rich cargo of proteins, RNAs, cytokines and lipids are selective ly packaged. Exosome could serve as an emerging platform for diagnostics. HQExo TM exosomes derived from body fluids could use as disease biomarkers research by ELISA, FACS, WB. Exosome are membranous struct ures and protected from degradation by extracellular proteases, which make the exosomes are highly stable. And due to the sample complexity reducing compared to the whole bodily fluids, exosome-based biomarker analy sis has attracted more and more attention to clinical diagnostics. Exosome can be purified by ultracentrifugation and characterized by nanoparticles tracking analysis (NTA) and ELISA or WB. Lyophilization is useful for a long-term storage at 4°C, and frozen liquid should be kept at -20°C to -80°C. Creative Biostructure standard exosome products guarantee higher purity and quality to meet our customer's downstream analyses.
Form	Lyophilized powder/ frozen liquid. Reconstitute lyophilized exosome by adding deionized water for a desired f inal concentration. Centrifuge before opening to ensure exosomes are at bottom, resuspend exosomes by pipett ing and/or vortex, please avoid bubbles. Centrifuge again and mix well for using.
Concentration	>1x10^8 particles
Storage	Lyophilized powder store at 4 °C. Frozen liquid store at -20°C to -80°C. Recommended to avoid repeated freez e-and-thaw cycles.