



PNExoTM Exosome-Magnolia

Catalog: PNE-FLM66

PRODUCT INFORMATION

Name	PNExo TM Exosome-Magnolia
Cat No.	PNE-FLM66
Source	Exosome derived from Magnolia
Product Overview	Plant exosomes are nanosized (30-150 nm) membrane vesicles that contain biomolecules. Plant-derived exosomes refer to naturally occurring nanoparticles derived from plants that contain bioactive molecules and proteins. These exosomes have been shown to have multiple benefits in a variety of applications, such as skincare, do godelivery, and biomedicine. Plant-derived exosomes have been found to possess antioxidant, anti-inflammated y, and anti-aging properties, making them an attractive option for the development of new and innovative their pies. Plant-derived natural substances are widely used as cosmeceutical materials because they exert beneficial effects on the human skin, such as antiaging, moisturizing, whitening, regeneration, and nutritional supply. But ides, they could delivery therapeutic compounds to target cells, potentially revolutionizing the way in which of ugs are administered. Overall, plant-derived exosomes hold great promise for a wide range of applications in effelds of medicine and biotechnology. PNExo TM is focused on the production and delivery of high quality plant-derived exosomes products. Exosomes are important tools of intercellular communication with a variety of biological functions, including cell regeneration and immune regulation. PNExo TM products undergo a rigoror s screening and purification process that guarantees their high purity and activity. Lyophilization is useful for long-term storage at 4°C, and frozen liquid should be kept at -20°C to -80°C. Ultracentrifugation and precipitation techniques are mainly used in exosome Isolation. It had been reported that both methods yielded extracell ar vesicles in the size range of exosomes and included apoproteins, which can be used in downstream analyse Creative Biostructure PNExo TM exosome products guarantee higher purity and quality to meet our customer research.
Form	Lyophilized powder
Concentration	> 1x10^6 particles
Storage	Lyophilized powder store at 4 °C. Frozen liquid store at -20°C to -80°C. Recommended to avoid repeated free e-and-thaw cycles.
Reconstitution	Reconstitute lyophilized exosome by adding deionized water for a desired final concentration. Centrifuge before opening to ensure exosomes are at bottom, resuspend exosomes by pipetting and/or vortex, please avoid but bles. Centrifuge again and mix well for using.

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel:1-631-317-1417 Fax:1-631-207-8356

Email: info@creative-biostructure.com