



Active uva Ursi Liposome

Catalog: Lipo-215C

This formulation contains liposomes made of phospholipids from soy lecithin and Uva ursi extract. The liposomal delivery system protects the active ingredient making it to reaches the target cells more efficiently. Uva ursi comes complete with other nutrients, including an
oxidants that provide natural protection against environmental stressors and have natural cleansing properties, which help minimize the appearance of pores, helping skin to appear smoother and more uniform while reducing puffiness. This product is intended to produce specific lotions creams and gels. Product designed to be used in cosmetic formulation in the range of 5 to 10 % respect to the final product.

PRODUCT INFORMATION

Product Name	Active uva Ursi Liposome
Catalog Number	Lipo-215C
Product Overview	This formulation contains liposomes made of phospholipids from soy lecithin and Uva ursi extract. The liposomal delivery system protects the active ingredient making it to reaches the target cells more efficiently. Uva ursi comes complete with other nutrients, including an oxidants that provide natural protection against environmental stressors and have natural cleansing properties, which help minimize the appearance of pores, helping skin to appear smoother and more uniform while reducing puffiness. This product is intended to produce specific lotions creams and gels. Product designed to be used in cosmetic formulation in the range of 5 to 10 % respect to the final product.
Active Ingredients	Uva Ursi Extract
Form	Viscous Nano Concentrates
Concentration	15%


Storage

Protected from light at room temperature.

USAGE GUIDELINES

General Handling	Avoid excessive mixing or vortexing to prevent aggregation. Long-term storage above -80°C may result in aggregate formation. Thaw on ice before use.
Storage	Short term: +2°C to +8°C (up to 1 week) Long term: -80°C (stable for 12 months)
Stability	Stable for 12 months at -80°C. Avoid repeated freeze-thaw cycles.
Freeze / Thaw	Can be frozen, but limit to 3 freeze/thaw cycles to maintain particle integrity.

Important: This product is for **research use only**. Not intended for diagnostic or therapeutic applications. Please refer to the Safety Data Sheet (SDS) for complete handling and safety information.