



Plain DPPC/CHOL/mPEG2000-DSPE liposomes (100 nm)

Catalog: CHOL-055L

This preformed liposome reagent is prepared in 10% sucrose and contains NO drug loading battery. It is typically used as a placebo control for the DPPC liposomes containing a drug loading battery, such as the DPPC liposomes containing ammonium sulfate gradient.
Lipid composition: DPPC/CHOL/mPEG2000-DSPE (50:45:5 mol/mol)

PRODUCT INFORMATION

Product Name	Plain DPPC/CHOL/mPEG2000-DSPE liposomes (100 nm)
Catalog Number	CHOL-055L
Product Overview	This preformed liposome reagent is prepared in 10% sucrose and contains NO drug loading battery. It is typically used as a placebo control for the DPPC liposomes containing a drug loading battery, such as the DPPC liposomes containing ammonium sulfate gradient. Lipid composition: DPPC/CHOL/mPEG2000-DSPE (50:45:5 mol/mol)
Lipid Composition	DPPC; CHOL; mPEG2000-DSPE
Application	Liposome production; Synthetic lipid
Storage Buffer	Internal/External solution: 10% sucrose, 20 mM HEPES, pH 7.0-7.5
Concentration	Lipid concentration: 60 mM (60 -65 mM)
Stability	6 Month for unopened vials.
Storage	Store at 2-8 centigrade.


Synonyms

DPPC; 16:0 PC,1,2-dipalmitoyl-sn-glycero-3-phosphocholine;
CHOL; cholesterol; mPEG2000-DSPE

USAGE GUIDELINES

General Handling	Avoid excessive mixing or vortexing to prevent aggregation. Long-term storage above -80°C may result in aggregate formation. Thaw on ice before use.
Storage	Short term: +2°C to +8°C (up to 1 week) Long term: -80°C (stable for 12 months)
Stability	Stable for 12 months at -80°C. Avoid repeated freeze-thaw cycles.
Freeze / Thaw	Can be frozen, but limit to ?3 freeze/thaw cycles to maintain particle integrity.

Important: This product is for **research use only**. Not intended for diagnostic or therapeutic applications. Please refer to the Safety Data Sheet (SDS) for complete handling and safety information.