



COVID-19 VLP (M; S; E Proteins)

Catalog: CBS-V952

PRODUCT INFORMATION

Cat No. CBS-V952 Source Insect cell recombinant Virus Name SARS-CoV-2 Virus Family Coronaviridae Composition Spike (S), envelope (E) and membrane (M) proteins Purity 85% Presentation 20mM Tris-Cl, pH7.9, 100mM NaCl and 0.5mM EDTA. Recombinant SARS-CoV-2 VLP that comprising Spike (S), envelope (E) and membrane (M) proteins was expression of the second		
Source Insect cell recombinant Virus Name SARS-CoV-2 Virus Family Coronaviridae Composition Spike (S), envelope (E) and membrane (M) proteins Purity 85% Presentation 20mM Tris-Cl, pH7.9, 100mM NaCl and 0.5mM EDTA. Recombinant SARS-CoV-2 VLP that comprising Spike (S), envelope (E) and membrane (M) proteins was expressed from insect cells. The VLPs were purified by a combination of gradient centrifugation and chromatograp	Name	COVID-19 VLP (M; S; E Proteins)
Virus Name SARS-CoV-2 Virus Family Coronaviridae Composition Spike (S), envelope (E) and membrane (M) proteins Purity 85% Presentation 20mM Tris-Cl, pH7.9, 100mM NaCl and 0.5mM EDTA. Recombinant SARS-CoV-2 VLP that comprising Spike (S), envelope (E) and membrane (M) proteins was expressed from insect cells. The VLPs were purified by a combination of gradient centrifugation and chromatographics.	Cat No.	CBS-V952
Virus Family Composition Spike (S), envelope (E) and membrane (M) proteins Purity 85% Presentation 20mM Tris-Cl, pH7.9, 100mM NaCl and 0.5mM EDTA. Recombinant SARS-CoV-2 VLP that comprising Spike (S), envelope (E) and membrane (M) proteins was expressed from insect cells. The VLPs were purified by a combination of gradient centrifugation and chromatographs.	Source	Insect cell recombinant
Composition Spike (S), envelope (E) and membrane (M) proteins Purity 85% Presentation 20mM Tris-Cl, pH7.9, 100mM NaCl and 0.5mM EDTA. Recombinant SARS-CoV-2 VLP that comprising Spike (S), envelope (E) and membrane (M) proteins was expensive to the second of the sec	Virus Name	SARS-CoV-2
Purity 85% Presentation 20mM Tris-Cl, pH7.9, 100mM NaCl and 0.5mM EDTA. Recombinant SARS-CoV-2 VLP that comprising Spike (S), envelope (E) and membrane (M) proteins was expection ressed from insect cells. The VLPs were purified by a combination of gradient centrifugation and chromatographs.	Virus Family	Coronaviridae
Presentation 20mM Tris-Cl, pH7.9, 100mM NaCl and 0.5mM EDTA. Recombinant SARS-CoV-2 VLP that comprising Spike (S), envelope (E) and membrane (M) proteins was exp Description ressed from insect cells. The VLPs were purified by a combination of gradient centrifugation and chromatographs.	Composition	Spike (S), envelope (E) and membrane (M) proteins
Recombinant SARS-CoV-2 VLP that comprising Spike (S), envelope (E) and membrane (M) proteins was exp Description ressed from insect cells. The VLPs were purified by a combination of gradient centrifugation and chromatograph	Purity	85%
Description ressed from insect cells. The VLPs were purified by a combination of gradient centrifugation and chromatograph	Presentation	20mM Tris-Cl, pH7.9, 100mM NaCl and 0.5mM EDTA.
	Description	Recombinant SARS-CoV-2 VLP that comprising Spike (S), envelope (E) and membrane (M) proteins was exp ressed from insect cells. The VLPs were purified by a combination of gradient centrifugation and chromatograp hic methods.

USAGE GUIDELINES

General	Avoid excessive mixing or shocking to prevent aggregation. Long term storage above -80°C may result in aggregate formation.
Storage	Short term: +2°C to +8°C Long term: -80°C
Stability	n.a.
Freezing	Can be frozen, but avoid multiple freeze/thaw cycles.

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel:1-631-317-1417 Fax:1-631-207-8356