



$HQExo^{TM}$ Exosome-SDH-Benign Prostatic Hyperplasia plasma

Catalog: Exo-HDBF-06

PRODUCT INFORMATION

Name	HQExo™ Exosome-SDH-Benign Prostatic Hyperplasia plasma
Cat No.	Exo-HDBF-06
Source	Exosome derived from Single Donor Human Benign Prostatic Hyperplasia plasma
Product Overview	Exosomes are nanosized vesicles (30-160 nm) secreted by exocytosis by most cell types and contain specifical cargos, such as RNAs, lipids, and proteins. The cargos amount and composition of exosomes depend on the collaboration of the collaboration of exosomes depend on the collaboration. Exosomes can deliver a variety of specific proteins, lipids and nucleic acids contained in them to near the protein of multiple physiological and pathological processes in the human body. Studies have shown that exosomes are related to the transport and release of characteristic molecules related to various diseases. The study of exosome from human disease-state body fluids will help us to systematically understand the relationship between exosomes and the occurrence and development of diseases. HQExo TM standard exosomes could use as positive controls for exosome isolation and functional research, such as ELISA, FACS, WB. Lyophilization is useful for a long-term storage at 4°C, and frozen liquid should be kept at -20°C to -80°C. Ultracentrifugation and precipitation techniques are mainly used in exosome Isolation. It had been reported that both methods yield ed extracellular vesicles in the size range of exosomes and included apoproteins, which can be used in downst eam analyses. Nanoparticles Tracking Analysis (NTA) is used for measuring exosome particles concentration, and WB or ELISA can be used in exosomal biomarkers analysis. Creative Biostructure standard exosome products guarantee higher purity and quality to meet our customer research.
Form	Lyophilized powder/ frozen liquid
Concentration	1x10^9 particles

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel:1-631-317-1417 Fax:1-631-207-8356