



DOPC-Based Fluorescent Liposome, Dextran-Rhodami ne B Encapsulated

Catalog: Lipo-051RG

PRODUCT INFORMATION

Name	DOPC-Based Fluorescent Liposome, Dextran-Rhodamine B Encapsulated
Cat No.	Lipo-051RG
Product Overview	As a large molecule, dextran can not leak out of the liposomes over time while the liposomes are inside the via l. By conjuncted with fluorescent probes, liposomes encapsulating fluorescent dextran can be used in experime nts that require tracking of the fluorescent water-soluble cargo. Excitation/Emission (nm): 570/590
Lipid Composition	DOPC/CHOL(7.0/3.0 mM) DOPC: 1,2-dioleoyl-sn-glycero-3-phosphocholine/CAS: 4235-95-4 CHOL: Choles terol/CAS: 57-88-5
Form	Liquid
Storage Buffer	20 mM HEPES, pH 7.4
Concentration	Lipid Concentration 10mM
Stability	2 months
Storage	2 °C to 8 °C (do not freeze)

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel:1-631-317-1417 Fax:1-631-207-8356