



POPE/CHOL/POPG Liposomes, Carboxyfluorescein E ncapsulated

Catalog: Lipo-058RG

PRODUCT INFORMATION

Name	POPE/CHOL/POPG Liposomes, Carboxyfluorescein Encapsulated
Cat No.	Lipo-058RG
Product Overview	Carboxyfluorescein is a self-quench dye and liposomes which encapsulate carboxyfluorescein can be used for full luorescent leakage assay. This liposome is made to order and delivered with no dialyzed and cleaned up. Custo mers should prepare the proper spin column for cleaning up the prep. It is advised to use the products within 30 days of the manufacturing date. Excitation/Emission (nm): 492/517 nm
Lipid Composition	POPE/CHOL/POPG (3.34/1.34/0.33 mM) POPE: 1-palmitoyl-2-oleoyl-sn-glycero-3-phosphoethanolamine/C AS:268550-95-4 CHOL: Cholesterol/CAS: 57-88-5 POPG: 1-palmitoyl-2-oleoyl-sn-glycero-3-phospho-(1'-ra c-glycerol) (sodium salt)/CAS: 268550-95-4
Form	Liquid
Storage Buffer	20 mM HEPES, pH 7.4
Concentration	Lipid Concentration 5mM
Stability	1 months
Storage	2 °C to 8 °C (do not freeze)

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel:1-631-317-1417 Fax:1-631-207-8356