



## **Mannosylated DiR Labeled Fluorescent Liposome, Con** trol

Catalog: Lipo-030RG

## PRODUCT INFORMATION

Name	Mannosylated DiR Labeled Fluorescent Liposome, Control
Cat No.	Lipo-030RG
Product Overview	DiR Labeled Fluorescent Liposome incorporate a lipophilic dye (DiR) inside their membranes. They are insoluble in water and the fluorescence is easily detected when incorporated into membranes. This liposome is mann osylated by 4-aminophenyl ?-D-mannopyranoside, which is a hydrophobic derivative of mannose. By mannosylation, the liposome could be targeted to the mannose recepter such as bacterial and viral cell surfaces. Excitation/Emission (nm): 750/780
Lipid Composition	L-alpha-PC/CHOL(24.3/10.9, mol/mol)
Form	Liquid
Storage Buffer	PBS, pH 7.4
Concentration	Lipid Concentration 35.1mM
Stability	2 months
Storage	2 °C to 8 °C (do not freeze)