



# DOTAP/DOPE(3/2 mM) Liposomes, Carboxyfluorescein Encapsulated

**Catalog: Lipo-090RG**

## PRODUCT INFORMATION

---

<b>Name</b>	DOTAP/DOPE(3/2 mM) Liposomes, Carboxyfluorescein Encapsulated
<b>Cat No.</b>	Lipo-090RG
<b>Product Overview</b>	Carboxy fluorescein is a self-quench dye and liposomes which encapsulate carboxyfluorescein can be used for fluorescent leakage assay. This liposome is made to order and delivered with no dialyzed and cleaned up. Customers should prepare the proper spin column for cleaning up the prep. It is advised to use the products within 30 days of the manufacturing date. Excitation/Emission (nm): 492/528 nm
<b>Lipid Composition</b>	DOTAP/DOPE(60/40 molar ratio)
<b>Form</b>	Liquid
<b>Storage Buffer</b>	20 mM HEPES, pH 7.4
<b>Concentration</b>	Lipid Concentration 5mM
<b>Stability</b>	1 months
<b>Storage</b>	2 °C to 8 °C (do not freeze)