



## DOTAP/DOPE(3/2 mM) Liposomes, Carboxyfluorescei n Encapsulated

Catalog: Lipo-090RG

## PRODUCT INFORMATION

Name	DOTAP/DOPE(3/2 mM) Liposomes, Carboxyfluorescein Encapsulated
Cat No.	Lipo-090RG
Product Overview	Carboxy fluorescein is a self-quench dye and liposomes which encapsulate carboxyfluorescein can be used for fluorescent leakage assay. This liposome is made to order and delivered with no dialyzed and cleaned up. Cust omers should prepare the proper spin column for cleaning up the prep. It is advised to use the products within 3 0 days of the manufacturing date. Excitation/Emission (nm): 492/528 nm
Lipid Composition	DOTAP/DOPE(60/40 molar ratio)
Form	Liquid
Storage Buffer	20 mM HEPES, pH 7.4
Concentration	Lipid Concentration 5mM
Stability	1 months
Storage	2 °C to 8 °C (do not freeze)

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel:1-631-317-1417 Fax:1-631-207-8356