



Negatively PS-Based Fluorescent Liposome, Dansyl lipi d (headgroup) Labeled

Catalog: Lipo-126RG

PRODUCT INFORMATION

Name	Negatively PS-Based Fluorescent Liposome, Dansyl lipid (headgroup) Labeled
Cat No.	Lipo-126RG
Product Overview	18:1 Dansyl PE has the Fluorescent group at the head group of the lipid. This product could be used for trackin g and detection. Excitation/Emission (nm): 336/513
Lipid Composition	L-alpha-PS/L-alpha-PC/CHOL(2/3/5 mM) L-alpha-PS: L-alpha-phosphatidylserine/CAS: 383908-63-2 L-alpha-PC: L-alpha-phosphatidylcholine/CAS: 97281-47-5 CHOL: Cholesterol/CAS: 57-88-5
Form	Liquid
Storage Buffer	PBS, pH 7.4
Concentration	Lipid Concentration 10mM
Stability	2 months
Storage	2 °C to 8 °C (do not freeze)