



# Lyophilized DOTAP/CHOL/DSPE-mPEG2000 ATP Liposome, Rhod-lipid Labeled

Catalog: Lipo-242RG

## PRODUCT INFORMATION

---

<b>Name</b>	Lyophilized DOTAP/CHOL/DSPE-mPEG2000 ATP Liposome, Rhod-lipid Labeled
<b>Cat No.</b>	Lipo-242RG
<b>Product Overview</b>	<p>The encapsulation of ATP in liposomes markedly promotes its effectiveness by preventing the hydrolysis by extracellular enzymes, increasing ATP circulation time and enhancing its intracellular penetration. ATP liposomes can be used in various models such as myocardial, liver, retina and wound healing ischemia. Studies have shown the ability of liposomal encapsulated ATP to prevent cell death and tissue dysfunction following ischemic events. The concentration of encapsulated ATP is 0.5<math>\mu</math>mol/vial.</p>
<b>Lipid Composition</b>	<p>DOTAP/CHOL/DSPE-mPEG2000/Rhod-PE (0.975/0.435/0.075/0.015 <math>\mu</math>mol/vial) DOTAP: 1,2-dioleoyl-3-trimethylammonium-propane (chloride salt) CHOL: Cholesterol DSPE-mPEG2000: 1,2-distearoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethyleneglycol)-2000] (ammonium salt)/CAS: 474922-77-5 Rhod PE: 1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-(lissamine rhodamine B sulfonyl) (ammonium salt) (Rhod PE)</p>
<b>Form</b>	Lyophilized Powder
<b>Storage Buffer</b>	PBS, pH 7.4 with trehalose as lyoprotectant
<b>Concentration</b>	Lipid Concentration 1.5 $\mu$ mol/vial
<b>Stability</b>	6 months
<b>Storage</b>	-20°C