



## DSPC/CHOL (55:45 mol/mol) Liposomes with Ammoni um Sulfate Gradient (100nm)

Catalog: DSPC-058L

## PRODUCT INFORMATION

Name	DSPC/CHOL (55:45 mol/mol) Liposomes with Ammonium Sulfate Gradient (100nm)
Cat No.	DSPC-058L
Product Overview	This product is prepared to the finest quality and shipped to you ready for drug loading. DSPC/CHO liposomes
	with ammonium gradient is the formulator's choice for loading of weakly basic compounds, because it offers si
	gnificantly improved drug retention and the ease of use in terms of process parameters such as loading tempera
	ture (60-70 deg C). One can use this liposomes as the placebo control for the drug loaded formulation or use th
	e much cheaper plain DSPC liposomes as placebo control.
	Lipid composition: DSPC/CHOL (55:45 mol/mol)
	Mean particle size: 100 nm (95 - 120 nm)
Lipid Composition	DSPC; CHOL
Application	Liposome production; Synthetic lipid
Storage Buffer	Hydration buffer (loading battery): 250 mM ammonium sulfate
	External solution: 10% sucrose, 20mM NaPO4, pH $6.5 \pm 0.2$
Concentration	Lipid concentration: 60 mM (60 - 70 mM)
Stability	6 Month for unopened vials.
Storage	Store at 2-8 centigrade.
Synonyms	DSPC; 18:0 PC,1,2-distearoyl-sn-glycero-3-phosphocholine; CHOL; cholesterol