



# Plain HSPC/CHOL liposomes (55:45 mol/mol, 100 nm)

**Catalog: CHOL-050L**

## PRODUCT INFORMATION

---

<b>Name</b>	Plain HSPC/CHOL liposomes (55:45 mol/mol, 100 nm)
<b>Cat No.</b>	CHOL-050L
<b>Product Overview</b>	<p>This preformed liposome reagent contains NO drug loading battery. It is typically used as a placebo control for the drug loaded liposomes with the same lipid composition.</p> <p>Lipid composition: HSPC/CHOL (55:45 mol/mol)</p> <p>Mean particle size: 100 nm (95 - 120 nm)</p>
<b>Lipid Composition</b>	HSPC; CHOL
<b>Application</b>	Liposome production; Synthetic lipid
<b>Storage Buffer</b>	Internal/External solution: 10% sucrose, 20 mM HEPES, pH 7.0 - 7.5
<b>Concentration</b>	Lipid concentration: 60 mM (60 -65 mM)
<b>Stability</b>	6 Month for unopened vials.
<b>Storage</b>	Store at 2-8 centigrade.
<b>Synonyms</b>	HSPC; hydrogenated phosphatidylcholine (Soy); CHOL; cholesterol