



## **EasyLoadTM LNP Encapsulation Kit (Lipid PreMix I)**

Catalog: Kit-BSLNP-001

## PRODUCT INFORMATION

Name	EasyLoadTM LNP Encapsulation Kit (Lipid PreMix I)
Cat No.	Kit-BSLNP-001
Product Overview	Creative Biostructure provides premix lipid kit for LNP encapsulation. The kit consists of ionizable lipids, base lipids, long-circulating lipids, and cholesterol. It can be used to prepare lipid nanoparticles loaded with nucleic acids, proteins, and drugs. The prepared lipid nanoparticles have good biocompatibility and high encapsulation efficiency, and are suitable for cell transfection and experimental animal injection. This kit contains one tube of lipid PreMix and one tube of citrate buffer. For Lipid PreMix I, the lipid components are: SM102, Cholestero l, DSPC, and PEG-2000-DMG. Please follow the instructions for the encapsulation process of LNP.
Composition	Lipid PreMix: 1mL Citrate buffer: 1mL
Form	Liquid
Concentration	Lipid PreMix: 16mM Citrate buffer: 50mM
Storage	Store at 4°C

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel:1-631-317-1417 Fax:1-631-207-8356

Email: info@creative-biostructure.com