



Placebo Control DOPC/CHOL Liposomes (100nm)

Catalog: DOPC-045L

This placebo control liposomes is made to match the lipid composition and characteristics for the non-PEGylated fluorescent liposomes. The product is sterile filtered and filled in sterile vials guaranteed to be stable for 6 months when stored in refrigerator. Lipid composition: DOPC/CHOL (55:45 mol/mol) Mean particle size: 100 nm (90-120 nm)

PRODUCT INFORMATION

Product Name	Placebo Control DOPC/CHOL Liposomes (100nm)
Catalog Number	DOPC-045L
Product Overview	This placebo control liposomes is made to match the lipid composition and characteristics for the non-PEGylated fluorescent liposomes. The product is sterile filtered and filled in sterile vials guaranteed to be stable for 6 months when stored in refrigerator. Lipid composition: DOPC/CHOL (55:45 mol/mol) Mean particle size: 100 nm (90-120 nm)
Lipid Composition	DOPC; CHOL
Application	Liposome production; Synthetic lipid

USAGE GUIDELINES

General Handling	Avoid excessive mixing or vortexing to prevent aggregation. Long-term storage above -80°C may result in aggregate formation. Thaw on ice before use.
Storage	Short term: +2°C to +8°C (up to 1 week) Long term: -80°C (stable for 12 months)



Stability	Stable for 12 months at -80°C. Avoid repeated freeze-thaw cycles.
Freeze / Thaw	Can be frozen, but limit to ?3 freeze/thaw cycles to maintain particle integrity.

Important: This product is for **research use only**. Not intended for diagnostic or therapeutic applications. Please refer to the Safety Data Sheet (SDS) for complete handling and safety information.