



Corning Plate 3552

Catalog: CB-111

PRODUCT INFORMATION

Product Name	Corning Plate 3552
Catalog Number	CB-111
Reservoir Well Volume	100 μ
Protein Well Reservoir	3
Protein Well Volume	2 μ
Material	Pzero
Surface	Pzero

USAGE GUIDELINES

General Handling	Avoid excessive mixing or vortexing to prevent aggregation. Long-term storage above -80°C may result in aggregate formation. Thaw on ice before use.
Storage	Short term: +2°C to +8°C (up to 1 week) Long term: -80°C (stable for 12 months)
Stability	Stable for 12 months at -80°C. Avoid repeated freeze-thaw cycles.
Freeze / Thaw	Can be frozen, but limit to $\leq$ 3 freeze/thaw cycles to maintain particle integrity.

Important: This product is for **research use only**. Not intended for diagnostic or therapeutic applications. Please refer to the Safety Data Sheet (SDS) for complete handling and safety information.